

Phospho-EIF2AK3 (Thr981) Antibody



PACO24119

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:**Gene ID:**

EIF2AK3

Uniprot

Q9NZJ5

Synonyms:

E2AK3; PERK; PEK; EIF2AK3;


Immunogen:

Peptide sequence around phosphorylation site of threonine 981(A-R-H(p)-T-G) derived from Human PEK/PERK.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded human prostate carcinoma tissue using PEK/PERK (Phospho-Thr982) antibody (left) or the same antibody preincubated with blocking peptide (right).