

# Phospho-RUNX1 (Ser435) Antibody



PACO24140

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:1000

**Protein Background:****Gene ID:**

RUNX1

**Uniprot**

Q01196

**Synonyms:**

Acute myeloid leukemia 1 protein; CBFA2; CBFA2; alpha 2 subunit; RUNX1

**Immunogen:**

Peptide sequence around phosphorylation site of Serine 435( S-N-S(p)-P-T) derived from Human AML1.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

N/A

N/A