

# Phospho-HCK (Tyr521) Antibody



PACO24165

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:1000

**Protein Background:****Gene ID:**

HCK

**Uniprot**

P08631

**Synonyms:**

p59-HCK/p60-HCK; Hemopoietic cell kinas; Tyrosine-protein kinase HCK;

**Immunogen:**


Peptide sequence around phosphorylation site of tyrosine 521(E-S-Q(p)-Y-Q) derived from Human HCK.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from JurKat cells using HCK (Phospho-Tyr522) Antibody. The lane on the right is treated with the antigen-specific peptide.