

Phospho-PPP1R2 (Ser120/Ser121) Antibody



PACO24212

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Inhibitor of protein-phosphatase 1.

Gene ID:

PPP1R2

Uniprot

P41236

Synonyms:

I-2; IPP-2; IPP2; Protein phosphatase inhibitor 2


Immunogen:

Peptide sequence around phosphorylation site of serine 120 and serine 121 (Q-E-S(p)-S(p)-GE) derived from Human PPP1R2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Immunohistochemistry analysis of paraffin-embedded human heart tissue using PPP1R2 (Phospho-Ser120+Ser121) antibody. The picture on the right is treated with the synthesized peptide.