

ALB Antibody



PACO24644

Product Information

Size:

50ug

Reactivity:

Sheep

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ALB

Uniprot

P14639

Synonyms:

Serum albumin, ALB

Immunogen:

Serum Albumin purified from Sheep serum (Native Protein).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A