

AADAC Antibody, HRP conjugated



PACO24673

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Displays cellular triglyceride lipase activity in liver, increases the levels of intracellular fatty acid, derived from the hydrolysis of newly formed triglyceride stores and plays a role in very low-density lipoprotein assembly. Displays serine esterase activity in liver. Deacetylates a variety of arylacetamide substrates, including xenobiotic compounds and procarcinogens, converting them to the primary arylamide compounds and increasing their toxicity.

Gene ID:

AADAC

Uniprot

P22760

Synonyms:

Arylacetamide deacetylase (EC 3.1.1.3), AADAC, DAC

Immunogen:

Recombinant Human Arylacetamide deacetylase protein (24-399AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A