

AADACL2 Antibody, Biotin conjugated



PACO24679

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	AADACL2
Source:	Uniprot
Rabbit	Q6P093
Isotype:	Synonyms:
IgG	Arylacetamide deacetylase-like 2 (EC 3.1.1. -), AADACL2
Applications:	Immunogen:
ELISA	Recombinant Human Arylacetamide deacetylase-like 2 protein (203-401AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A