

ACADL Antibody, Biotin conjugated



PACO24707

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

mitochondrial matrix, acyl-CoA dehydrogenase activity, electron carrier activity, fatty-acyl-CoA binding, flavin adenine dinucleotide binding, long-chain-acyl-CoA dehydrogenase activity, oxidoreductase activity, acting on the CH-CH group of donors, with a flavin as acceptor, palmitoyl-CoA oxidase activity, carnitine catabolic process, carnitine metabolic process, CoA-linked.

Gene ID:

ACADL

Uniprot

P28330

Synonyms:

Long-chain specific acyl-CoA dehydrogenase, mitochondrial (LCAD) (EC 1.3.8.8), ACADL

Immunogen:

Recombinant Human Long-chain specific acyl-CoA dehydrogenase, mitochondrial protein (30-260AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A