

ACRV1 Antibody



PACO24728

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

multicellular organismal development.

Gene ID:

ACRV1

Uniprot

P26436

Synonyms:

Acrosomal protein SP-10 (Acrosomal vesicle protein 1), ACRV1


Immunogen:

Recombinant Human Acrosomal protein SP-10 protein (22-265AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO24728 at dilution of 1:100.