

ACRV1 Antibody, Biotin conjugated



PACO24731

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

multicellular organismal development.

Gene ID:

ACRV1

Uniprot

P26436

Synonyms:

Acrosomal protein SP-10 (Acrosomal vesicle protein 1), ACRV1

Immunogen:

Recombinant Human Acrosomal protein SP-10 protein (22-265AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A