

Afp Antibody, Biotin conjugated



PACO24802

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds estrogens, fatty acid, and metals.

Gene ID:

Afp

Uniprot

P02772

Synonyms:

Alpha-fetoprotein (Alpha-1-fetoprotein) (Alpha-fetoglobulin), Afp

Immunogen:

Recombinant Mouse Alpha-fetoprotein protein (18-605AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A