

# AFTPH Antibody



PACO24803

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May play a role in membrane trafficking.

**Gene ID:**

AFTPH

**Uniprot**

Q6ULP2

**Synonyms:**

Aftiphilin, AFTPH, AFTH

**Immunogen:**


Recombinant Human Aftiphilin protein (752-889AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human brain tissue using PACO24803 at dilution of 1:100.