

AGA Antibody, Biotin conjugated



PACO24810

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cleaves the GlcNAc-Asn bond which joins oligosaccharides to the peptide of asparagine-linked glycoproteins.

Gene ID:

AGA

Uniprot

P20933

Synonyms:

N(4)-(beta-N-acetylglucosaminy)-L-asparaginase (EC 3.5.1.26)
(Aspartylglucosaminidase) (Glycosylasparaginase) (N4-(N-acetyl-beta-glucosaminy)-L-asparagine amidase) [Cleaved into: Glycosylasparaginase alpha chain; Glycosylasparaginase beta chain], AGA

Immunogen:

Recombinant Human N(4)-(beta-N-acetylglucosaminy)-L-asparaginase protein (206-346AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A