

PACO24811

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Dual metabolic roles of gluconeogenesis (in the mitochondria) and glyoxylate detoxification (in the peroxisomes).

Gene ID:

Agxt

Uniprot

O35423

Synonyms:

Serine--pyruvate aminotransferase, mitochondrial (SPT) (EC 2.6.1.51) (Alanine--glyoxylate aminotransferase) (AGT) (EC 2.6.1.44), Agxt, Agxt1

Immunogen:

Recombinant Mouse Serine--pyruvate aminotransferase, mitochondrial protein (25-414AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A