

# AK1 Antibody, Biotin conjugated



PACO24838

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Catalyzes the reversible transfer of the terminal phosphate group between ATP and AMP. Also displays broad nucleoside diphosphate kinase activity. Plays an important role in cellular energy homeostasis and in adenine nucleotide metabolism.

**Gene ID:**

AK1

**Uniprot**

P00568

**Synonyms:**

Adenylate kinase isoenzyme 1 (AK 1) (EC 2.7.4.3) (EC 2.7.4.6) (ATP-AMP transphosphorylase 1) (ATP: AMP phosphotransferase) (Adenylate monophosphate kinase) (Myokinase), AK1

**Immunogen:**

Recombinant Human Adenylate kinase isoenzyme 1 protein (1-194AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A