

ANXA9 Antibody, FITC conjugated



PACO24942

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Low affinity receptor for acetylcholine known to be targeted by disease-causing pemphigus vulgaris antibodies in keratinocytes.

Gene ID:

ANXA9

Uniprot

O76027

Synonyms:

Annexin A9 (Annexin XXXI) (Annexin-31) (Annexin-9) (Pemphaxin), ANXA9, ANX31

Immunogen:

Recombinant Human Annexin A9 protein (8-345AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A