

# AP1S1 Antibody, FITC conjugated



PACO24950

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Subunit of clathrin-associated adaptor protein complex 1 that plays a role in protein sorting in the late-Golgi/trans-Golgi network (TGN) and/or endosomes. The AP complexes mediate both the recruitment of clathrin to membranes and the recognition of sorting signals within the cytosolic tails of transmembrane cargo molecules.

**Gene ID:**

AP1S1

**Uniprot**

P61966

**Synonyms:**

AP-1 complex subunit sigma-1A (Adaptor protein complex AP-1 subunit sigma-1A) (Adaptor-related protein complex 1 subunit sigma-1A) (Clathrin assembly protein complex 1 sigma-1A small chain) (Clathrin coat assembly protein AP19) (Golgi adaptor HA1/AP1 adaptin sigma-1A subunit) (HA1 19 kDa subunit) (Sigma 1a subunit of AP-1 clathrin) (Sigma-adaptin 1A) (Sigma1A-adaptin), AP1S1, AP19 CLAPS1

**Immunogen:**

Recombinant Human AP-1 complex subunit sigma-1A protein (1-127AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A