

AP3S2 Antibody, HRP conjugated



PACO24961

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Part of the AP-3 complex, an adaptor-related complex which is not clathrin-associated. The complex is associated with the Golgi region as well as more peripheral structures. It facilitates the budding of vesicles from the Golgi membrane and may be directly involved in trafficking to lysosomes. In concert with the BLOC-1 complex, AP-3 is required to target cargos into vesicles assembled at cell bodies for delivery into neurites and nerve terminals.

Gene ID:

AP3S2

Uniprot

P59780

Synonyms:

AP-3 complex subunit sigma-2 (AP-3 complex subunit sigma-3B) (Adaptor-related protein complex 3 subunit sigma-2) (Sigma-3B-adaptin) (Sigma3B-adaptin) (Sigma-adaptin 3b), AP3S2

Immunogen:

Recombinant Human AP-3 complex subunit sigma-2 protein (1-193AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A