

APAF1 Antibody, Biotin conjugated



PACO24967

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Oligomeric Apaf-1 mediates the cytochrome c-dependent autocatalytic activation of pro-caspase-9 (Apaf-3), leading to the activation of caspase-3 and apoptosis. This activation requires ATP. Isoform 6 is less effective in inducing apoptosis.

Gene ID:

APAF1

Uniprot

O14727

Synonyms:

Apoptotic protease-activating factor 1 (APAF-1), APAF1, KIAA0413

Immunogen:

Recombinant Human Apoptotic protease-activating factor 1 protein (104-415AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A