

Apoa4 Antibody, HRP conjugated



PACO24991

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May have a role in chylomicrons and VLDL secretion and catabolism. Required for efficient activation of lipoprotein lipase by ApoC-II; potent activator of LCAT. ApoA-IV is a major component of HDL and chylomicrons.

Gene ID:

Apoa4

Uniprot

P06728

Synonyms:

Apolipoprotein A-IV (Apo-AIV) (ApoA-IV) (Apolipoprotein A4), Apoa4

Immunogen:

Recombinant Mouse Apolipoprotein A-IV protein (21-395AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A