

# Apoa4 Antibody, FITC conjugated



PACO24992

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May have a role in chylomicrons and VLDL secretion and catabolism. Required for efficient activation of lipoprotein lipase by ApoC-II; potent activator of LCAT. ApoA-IV is a major component of HDL and chylomicrons.

**Gene ID:**

Apoa4

**Uniprot**

P06728

**Synonyms:**

Apolipoprotein A-IV (Apo-AIV) (ApoA-IV) (Apolipoprotein A4), Apoa4

**Immunogen:**

Recombinant Mouse Apolipoprotein A-IV protein (21-395AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A