

APOA4 Antibody, HRP conjugated



PACO24995

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May have a role in chylomicrons and VLDL secretion and catabolism. Required for efficient activation of lipoprotein lipase by ApoC-II; potent activator of LCAT. ApoA-IV is a major component of HDL and chylomicrons.

Gene ID:

APOA4

Uniprot

P06727

Synonyms:

Apolipoprotein A-IV (Apo-AIV) (ApoA-IV) (Apolipoprotein A4), APOA4

Immunogen:

Recombinant Human Apolipoprotein A-IV protein (21-396AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A