

APOD Antibody, FITC conjugated



PACO25008

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

APOD occurs in the macromolecular complex with lecithin-cholesterol acyltransferase. It is probably involved in the transport and binding of bilin. Appears to be able to transport a variety of ligands in a number of different contexts.

Gene ID:

APOD

Uniprot

P05090

Synonyms:

Apolipoprotein D (Apo-D) (ApoD), APOD

Immunogen:

Recombinant Human Apolipoprotein D protein (21-189AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A