

ApoE Antibody, Biotin conjugated



PACO25013

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Mediates the binding, internalization, and catabolism of lipoprotein particles. It can serve as a ligand for the LDL (apo B/E) receptor and for the specific apo-E receptor (chylomicron remnant) of hepatic tissues.

Gene ID:

ApoE

Uniprot

P08226

Synonyms:

Apolipoprotein E (Apo-E), ApoE

Immunogen:

Recombinant Mouse Apolipoprotein E protein (19-311AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A