

ARFIP2 Antibody, Biotin conjugated



PACO25027

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Putative target protein of ADP-ribosylation factor. Involved in membrane ruffling.

Gene ID:

ARFIP2

Uniprot

P53365

Synonyms:

Arfaptin-2 (ADP-ribosylation factor-interacting protein 2) (Partner of RAC1) (Protein POR1), ARFIP2, POR1

Immunogen:

Recombinant Human Arfaptin-2 protein (20-295AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A