

ARL4D Antibody, Biotin conjugated



PACO25051

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Small GTP-binding protein which cycles between an inactive GDP-bound and an active GTP-bound form, and the rate of cycling is regulated by guanine nucleotide exchange factors (GEF) and GTPase-activating proteins (GAP). GTP-binding protein that does not act as an allosteric activator of the cholera toxin catalytic subunit. Recruits CYTH1, CYTH2, CYTH3 and CYTH4 to the plasma membrane in GDP-bound form.

Gene ID:

ARL4D

Uniprot

P49703

Synonyms:

ADP-ribosylation factor-like protein 4D (ADP-ribosylation factor-like protein 4L), ARL4D, ARF4L

Immunogen:

Recombinant Human ADP-ribosylation factor-like protein 4D protein (2-201AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A