

PACO25064

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:**Gene ID:**

ARSB

Uniprot

P15848

Synonyms:

Arylsulfatase B (ASB) (EC 3.1.6.12) (N-acetylgalactosamine-4-sulfatase) (G4S), ARSB


Immunogen:

Recombinant Human Arylsulfatase B protein (37-533AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ARSB antibody at 2 μ g/ml + Mouse brain tissue

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 60, 46 kDa

Observed band size: 40 kDa