

ASL Antibody, HRP conjugated



PACO25093

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytoplasm, cytosol, extracellular exosome, argininosuccinate lyase activity, arginine biosynthetic process, arginine catabolic process, cellular nitrogen compound metabolic process, internal protein amino acid, acetylation, polyamine metabolic process, protein tetramerization.

Gene ID:

ASL

Uniprot

P04424

Synonyms:

Argininosuccinate lyase (ASAL) (EC 4.3.2.1) (Arginosuccinase), ASL

Immunogen:

Recombinant Human Argininosuccinate lyase protein (2-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A