

# ASL Antibody, Biotin conjugated



PACO25095

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytoplasm, cytosol, extracellular exosome, argininosuccinate lyase activity, arginine biosynthetic process, arginine catabolic process, cellular nitrogen compound metabolic process, internal protein amino acid, acetylation, polyamine metabolic process, protein tetramerization.

**Gene ID:**

ASL

**Uniprot**

P04424

**Synonyms:**

Argininosuccinate lyase (ASAL) (EC 4.3.2.1) (Arginosuccinase), ASL

**Immunogen:**

Recombinant Human Argininosuccinate lyase protein (2-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A