

# CSTB Antibody, FITC conjugated



PACO25102

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is a respiratory chain that generates an electrochemical potential coupled to ATP synthesis.

**Gene ID:**

CSTB

**Uniprot**

P04080

**Synonyms:**

Cystatin-B (CPI-B) (Liver thiol proteinase inhibitor) (Stefin-B), CSTB, CST6 STFB

**Immunogen:**

Recombinant Human Cystatin-B protein (1-98AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A