

CSTB Antibody, Biotin conjugated



PACO25103

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is a respiratory chain that generates an electrochemical potential coupled to ATP synthesis.

Gene ID:

CSTB

Uniprot

P04080

Synonyms:

Cystatin-B (CPI-B) (Liver thiol proteinase inhibitor) (Stefin-B), CSTB, CST6 STFB

Immunogen:

Recombinant Human Cystatin-B protein (1-98AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A