

ATP6V1G2 Antibody, HRP conjugated



PACO25157

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalytic subunit of the peripheral V1 complex of vacuolar ATPase (V-ATPase). V-ATPase is responsible for acidifying a variety of intracellular compartments in eukaryotic cells.

Gene ID:

ATP6V1G2

Uniprot

O95670

Synonyms:

V-type proton ATPase subunit G 2 (V-ATPase subunit G 2) (V-ATPase 13 kDa subunit 2) (Vacuolar proton pump subunit G 2), ATP6V1G2, ATP6G ATP6G2 NG38

Immunogen:

Recombinant Human V-type proton ATPase subunit G 2 protein (1-118AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A