

ATPAF2 Antibody, HRP conjugated



PACO25165

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	ATPAF2
Source:	Uniprot
Rabbit	Q8N5M1
Isotype:	Synonyms:
IgG	ATP synthase mitochondrial F1 complex assembly factor 2 (ATP12 homolog), ATPAF2, ATP12
Applications:	Immunogen:
ELISA	Recombinant Human ATP synthase mitochondrial F1 complex assembly factor 2 protein (41-289AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A