

P4HB Antibody, HRP conjugated



PACO25209

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Dsb proteins (DsbA, DsbB, DsbC, and DsbD) catalyze formation and isomerization of protein disulfide bonds in the periplasm of Escherichia coli. Disulphide bond isomerase (DsbC), a member of the thioredoxin superfamily, acts as a disulfide isomerase and converts aberrant disulfide bonds to correct ones. It also acts by transferring its disulfide bond to other proteins and is reduced in the process. DsbC is reoxidized by a yet uncharacterized protein.

Gene ID:

P4HB

Uniprot

P07237

Synonyms:

Protein disulfide-isomerase (PDI) (EC 5.3.4.1) (Cellular thyroid hormone-binding protein) (Prolyl 4-hydroxylase subunit beta) (p55), P4HB, ERBA2L PDI PDIA1 PO4DB

Immunogen:

Recombinant Human Protein disulfide-isomerase protein (29-508AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A