

BDH1 Antibody, FITC conjugated



PACO25246

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

D-beta-Hydroxybutyrate dehydrogenase mitochondria has been purified to apparent homogeneity. The membrane-bound enzyme is first released by phospholipase A digestion of the mitochondria.

Gene ID:

BDH1

Uniprot

Q02338

Synonyms:

D-beta-hydroxybutyrate dehydrogenase, mitochondrial (EC 1.1.1.30) (3-hydroxybutyrate dehydrogenase) (BDH) (Short chain dehydrogenase/reductase family 9C member 1), BDH1, BDH SDR9C1

Immunogen:

Recombinant Human D-beta-hydroxybutyrate dehydrogenase, mitochondrial protein (47-343AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A