

# Bglap Antibody, Biotin conjugated



PACO25255

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Constitutes 1-2% of the total bone protein. It binds strongly to apatite and calcium.

**Gene ID:**

Bglap

**Uniprot**

P04640

**Synonyms:**

Osteocalcin (Bone Gla protein) (BGP) (Gamma-carboxyglutamic acid, containing protein), Bglap, Bglap2

**Immunogen:**

Recombinant Rat Osteocalcin protein (50-99AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A