

BGN Antibody, Biotin conjugated



PACO25259

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in collagen fiber assembly.

Gene ID:

BGN

Uniprot

P21810

Synonyms:

Biglycan (Bone/cartilage proteoglycan I) (PG-S1), BGN, SLRR1A

Immunogen:

Recombinant Human Biglycan protein (38-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A