

BMP3 Antibody, HRP conjugated



PACO25269

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Could be involved in bone and cartilage formation. Chondrogenic signaling is mediated by the high-affinity receptor BMPR1B.

Gene ID:

BMP3

Uniprot

P12645

Synonyms:

Bone morphogenetic protein 3 (BMP-3) (Bone morphogenetic protein 3A) (BMP-3A) (Osteogenin), BMP3, BMP3A

Immunogen:

Recombinant Human Bone morphogenetic protein 3 protein (362-472AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A