

# BMP3 Antibody, FITC conjugated



PACO25270

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Could be involved in bone and cartilage formation. Chondrogenic signaling is mediated by the high-affinity receptor BMPR1B.

**Gene ID:**

BMP3

**Uniprot**

P12645

**Synonyms:**

Bone morphogenetic protein 3 (BMP-3) (Bone morphogenetic protein 3A) (BMP-3A) (Osteogenin), BMP3, BMP3A

**Immunogen:**

Recombinant Human Bone morphogenetic protein 3 protein (362-472AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A