

LGB Antibody, HRP conjugated



PACO25285

Product Information

Size:

50ug

Reactivity:

Bovine

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

LGB

Uniprot

P02754

Synonyms:

Beta-lactoglobulin (Beta-LG) (allergen Bos d 5), LGB

Immunogen:

Bovine milk Beta-lactoglobulin (Native Protein).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A