

BTN3A2 Antibody, FITC conjugated



PACO25294

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role in T-cell responses in the adaptive immune response. Inhibits the release of IFNG from activated T-cells.

Gene ID:

BTN3A2

Uniprot

P78410

Synonyms:

Butyrophilin subfamily 3 member A2, BTN3A2, BT3.2 BTF3 BTF4

Immunogen:

Recombinant Human Butyrophilin subfamily 3 member A2 protein (30-248AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A