

# BTN3A2 Antibody, Biotin conjugated



PACO25295

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Plays a role in T-cell responses in the adaptive immune response. Inhibits the release of IFNG from activated T-cells.

**Gene ID:**

BTN3A2

**Uniprot**

P78410

**Synonyms:**

Butyrophilin subfamily 3 member A2, BTN3A2, BT3.2 BTF3 BTF4

**Immunogen:**

Recombinant Human Butyrophilin subfamily 3 member A2 protein (30-248AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A