

BTN3A3 Antibody, Biotin conjugated



PACO25299

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role in T-cell responses in the adaptive immune response.

Gene ID:

BTN3A3

Uniprot

O00478

Synonyms:

Butyrophilin subfamily 3 member A3, BTN3A3, BTF3

Immunogen:

Recombinant Human Butyrophilin subfamily 3 member A3 protein (30-248AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A