

LINC01553 Antibody



PACO25304

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

LINC01553

Uniprot

A4QN01

Synonyms:

Putative uncharacterized protein encoded by LINC01553, LINC01553, C10orf40


Immunogen:

Recombinant Human Putative uncharacterized protein encoded by LINC01553 protein (1-128AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: LINC01553 antibody at 12 μ g/ml

Lane 1: Rat heart tissue

Lane 2: Mouse skeletal muscle tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 16 kDa

Observed band size: 16 kDa

.