

SMCO3 Antibody, HRP conjugated



PACO25313

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SMCO3

Uniprot

A2RU48

Synonyms:

Single-pass membrane and coiled-coil domain-containing protein 3, SMCO3, C12orf69

Immunogen:

Recombinant Human Single-pass membrane and coiled-coil domain-containing protein 3 protein (1-154AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A