

**alpha-1-acidglycoprotein Antibody, HRP
conjugated
PACO25317**



Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

Synonyms:

Immunogen:

Human α alpha-1-acidglycoprotein (Native Protein).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A