

# RAB1A Antibody, FITC conjugated



PACO25322

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Probably required for transit of protein from the ER through Golgi compartment. Binds GTP and GDP and possesses intrinsic GTPase activity.

**Gene ID:**

RAB1A

**Uniprot**

P62820

**Synonyms:**

Ras-related protein Rab-1A (YPT1-related protein), RAB1A, RAB1

**Immunogen:**

Recombinant Human Ras-related protein Rab-1A protein (2-205AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A