

RAB1A Antibody, Biotin conjugated



PACO25323

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probably required for transit of protein from the ER through Golgi compartment. Binds GTP and GDP and possesses intrinsic GTPase activity.

Gene ID:

RAB1A

Uniprot

P62820

Synonyms:

Ras-related protein Rab-1A (YPT1-related protein), RAB1A, RAB1

Immunogen:

Recombinant Human Ras-related protein Rab-1A protein (2-205AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A