

TTR Antibody, FITC conjugated



PACO25330

Product Information

Size:

50ug

Reactivity:

Human

Source:

Guinea pig

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain. Detected in serum and cerebrospinal fluid (at protein level). Highly expressed in choroid plexus epithelial cells. Detected in retina pigment epithelium and liver.

Gene ID:

TTR

Uniprot

P02766

Synonyms:

Transthyretin (ATTR) (Prealbumin) (TBPA), TTR, PALB

Immunogen:

Human Prealbumin protein (Native Protein).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A