

PACO25332

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

IMUP

Uniprot

Q9GZP8

Synonyms:

Immortalization up-regulated protein (Hepatocyte growth factor activator inhibitor type 2-related small protein) (H2RSP) (HAI-2-related small protein), IMUP, C19orf33

Immunogen:

Recombinant Human Immortalization up-regulated protein (1-85AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A